

HPLC Injector

Lab Resources

Posted by: webmaster

Posted on : 2012/8/17 18:42:17

HPLC Injector

Rheodyne



1. Manual Injector

Rheodyne HPLC Sample Manual Injector

Cat. No. Description

70102-position, 6-port, single mode, SS

77252-position, 6-port, dual mode with MBB, SS

7725i2-position, 6-port, dual mode with MBB and position sensor, SS

90102-position, 6-port, single mode, PK

97252-position, 6-port, dual mode with MBB and position sensor, PK

9725i2-position, 6-port, dual mode with MBB and position sensor, PK

2. RheBuild Kits

Cat. No. Description

7010-996RheBuild Conversion Kit to include Stator Face Assembly for 7010

7010-997RheBuild Kit w/Stator for 7010

7010-998RheBuild pH Upgrade Kit for 7000 Series

7010-999RheBuild Kit for 7010/7000

7125-998RheBuild pH Upgrade Kit for 7125

7125-999RheBuild Kit for 7125

7125Ti-999RheBuild Kit for 7125-081

7410-999RheBuild Kit for 7410 Rheodyne® Valves

7520-999RheBuild Kit for 7520

77104-999RheBuild Kit for 7710-040

7725-999RheBuild Kit for 7725i

7739-999RheBuild Kit for 7739

7741-999RheBuild Kit for 7710-040

8125-999RheBuild Kit for 8125/8126

9010-999RheBuild Kit for 9010/9040

9125-999RheBuild Kit for 9125/9126

9725-999RheBuild Kit for 9725/9725



3. Stainless Steel Sample Loops

Loops for 7125 and 7010 Injectors

Cat. No. Description

70205µL SS Sample Loop for 7010/7125

702110µL SS Sample Loop for 7010/7125

702220µL SS Sample Loop for 7010/7125

702350µL SS Loop for 7010/7125

7024100µL SS Loop for 7010/7125

7025200µL SS Loop for 7010/7125

7026500µL SS Sample Loop for 7010/7125

70271mL SS Loop for 7010/7125

70282mL SS Loop for 7010/7125

70295mL SS Sample Loop for 7010/7125



Loops for 7725(i), PR700-100, and PR703-100

Cat. No. Description

7755-0205µL SS Loop,
7755-02110µL SS Loop,
7755-02220µL SS Sample Loop,
7755-02350µL SS Sample Loop,
7755-024100µL SS Loop,
7755-025200µL SS Sample Loop,
7755-026500µL SS Sample Loop,
7755-0271mL SS Loop,
7755-0295mL SS Loop



Loops for 3725(i)-038 Injectors

Cat. No. Description

3065-0182mL Prep Sample Loop SS
3065-0195mL Prep Sample Loop SS
3065-02310mL Prep Sample Loop SS
3065-02520mL Prep Sample Loop SS

Loops For 8125 Injectors

Cat. No. Description

80205µL SS Sample Loop for 8125/8126
802110µL SS Loop for 8125/8126
802220µL SS Sample Loop for 8125/8126
802350µL SS Sample Loop for 8125/8126

4. Rotor Seals

Loops For 8125 Injectors

Cat. No. Description

80205µL SS Sample Loop for 8125/8126
802110µL SS Loop for 8125/8126
802220µL SS Sample Loop for 8125/8126
802350µL SS Sample Loop for 8125/8126



Tefzel Rotor Seals Tefzel blend rotor seals have an operating pH range from 0 to 14.

Cat. No. Description

7010-071 Tefzel Rotor Seal for 7010, 7010-087, 7000, 7040

7030-015Tefzel Rotor Seal for 7030 and 9030
7060-074Tefzel Rotor Seal for 7060, 7066, 9060
7125-079Tefzel Rotor Seal for 7125, 7125-081, 7725
7410-075Tefzel Rotor Seal for 7410
8125-097Tefzel Rotor Seal for 8125, 8126
9010-051Tefzel Rotor Seal for 9010
9125-082Tefzel Rotor Seal for 9125, 9725

Vespel Rotor SealsVespel blend rotor seals have an operating pH range from 0 to 10.

Cat. No.Description

7010-039Vespel Rotor Seal for 7000, 7010, 7040
7030-003Vespel Rotor Seal for 7030
7060-070Vespel Rotor Seal for 7060, 7066
7125-047Vespel Rotor Seal for 7125, 7725
7410-038Vespel Rotor Seal for 7410
7413-013Vespel Rotor Seal for 7413
8125-038Vespel Rotor Seal for 8125

5. Stators

DuraLife® and UltraLife(TM) Stators

Cat. No.Description

Stainless Steel Stators

3725-085SS Stator for 3725-038, 3710-038, 3000-038, 3030-038
7010-040SS Stator for 7010, 7125, 7000, 7030, 7040
7060-039SS Stator for 7060, 7066
7123-047SS Stator for EV/PR/RV500-100
7123-128SS Stator for EV/HT/PR/RV700-107
7123-142SS Stator for EV/PR/RV500-104
7123-148SS Stator for EV/PR/RV500-101
7123-180SS Stator for HT/PD703-100, EV/RV703-100
7123-223SS Stator for EV/HT/RV700-112
7410-041SS Stator for 7410, 7413
7650-002SS Stator for EV/PR/RV700-102
7725-010SS Stator for 7725
7750-038SS Stator for EV/PR/RV700-100
7900-107SS Stator for MX7900-000
7925-002SS Stator for MX7925-000
7960-002SS Stator for MX7960-000
8125-098SS Stator for 8125, 8126



